



June 3, 2016

To: Directors of Clinical Laboratories
Directors of Transfusion Services

Subject: RBC Components from Apheresis Collections

LifeSouth is implementing new technology and expanding its options for red blood cell collections by apheresis. The lists below include new leukoreduced and non-leukoreduced ISBT product codes.

Please note: your computer system may require updates in order to accept these new components.

Leukoreduced Product Codes	
E4543	Apheresis RED BLOOD CELLS ACD-A>AS3/XX/refg ResLeu:<5log6
E4538	Apheresis RED BLOOD CELLS ACD-A>AS3/XX/refg Irradiated ResLeu:<5log6
E4587	Deglycerolized Apheresis RED BLOOD CELLS None/XX/refg Open ResLeu:<5E6
E4582	Deglycerolized Apheresis RED BLOOD CELLS None/XX/refg Open Irradiated ResLeu:<5E6

Non-leukoreduced Product Codes	
E4524	Apheresis RED BLOOD CELLS ACD-A>AS1/XX/refg
E4534	Apheresis RED BLOOD CELLS ACD-A>AS1/XX/refg 1st container
E4535	Apheresis RED BLOOD CELLS ACD-A>AS1/XX/refg 2nd container
E4525	Apheresis RED BLOOD CELLS ACD-A>AS1/XX/refg Irradiated
E4529	Apheresis RED BLOOD CELLS ACD-A>AS1/XX/refg Irradiated 1st container
E4530	Apheresis RED BLOOD CELLS ACD-A>AS1/XX/refg Irradiated 2nd container
E4537	Apheresis RED BLOOD CELLS ACD-A>AS3/XX/refg Irradiated
E4580	Deglycerolized Apheresis RED BLOOD CELLS None/XX/refg Open
E4590	Deglycerolized Apheresis RED BLOOD CELLS None/XX/refg Open 1st container
E4591	Deglycerolized Apheresis RED BLOOD CELLS None/XX/refg Open 2nd container
E4581	Deglycerolized Apheresis RED BLOOD CELLS None/XX/refg Open Irradiated
E4585	Deglycerolized Apheresis RED BLOOD CELLS None/XX/refg Open Irradiated 1st container
E4586	Deglycerolized Apheresis RED BLOOD CELLS None/XX/refg Open Irradiated 2nd container

Please feel free to contact the blood center's medical director with any questions or concerns you have about these products: Juan A. Merayo-Rodriguez, MD, 1-888-795-2707, jamerayo-rodriguez@lifesouth.org